



## Centrifuge "Centro-4-BL"

WITH SWINGOUT HEAD FOR 4 TUBES OF 75 X 13 mm VAC.  
FREE MAINTENANCE INDUCTION DRIVE MOTOR.



Control panel features (see page 116).

Centrifuge kit for separating plasma rich platelets from whole blood.



### MODEL

Part No.	Height / Width / Depth (exterior) cm	Power W	Max. R.p.m.	Max. R.c.f. (xg)	Weight Kg
<b>7002560</b>	20 28 34	170	5000	2907	6

Rotor comes without tubes.

ABS plastic moulded external case.



## Centrifuge "Centro-8-BL"

WITH ANGLE ROTOR FOR 7-15 ml TUBES AND 75 X 13 mm VAC (WITH ADAPTERS).  
FREE MAINTENANCE INDUCTION DRIVE MOTOR.



Control panel features (see page 116).



ABS plastic moulded external case



### MODEL

Part No.	Max. Volume ml.	Tube capacity	Height / Width / Depth (exterior) cm	Power W	Max. R.p.m.	Max. R.c.f. (xg)	Weight Kg
<b>7002356</b>	120	8 x 15 ml	20 28 34	170	4500	2515	6.5

Supplied complete with 8 x 15 ml glass tubes.

### ACCESSORIES

**Bucket adapter for 7 ml tubes.**  
Rotor capacity: 8 adapters.  
Part No. **7000700**

**Bucket adapter for 75 x 13 mm VAC. tubes.**  
Rotor capacity: 8 adapters.  
Part No. **1001557**

**SPARES. Tubes** (see page 132).